


Assigning an active tag to an animal







1. Go to **> My Devices > Add device > Active Tag** to begin pairing.
2. When prompted, swipe the magnet over the magnet icon until the LED flashes TWICE. Follow the instructions to pair the animal VID/LID or EID.



Attaching a collar tag to an animal

1. Position weight and active tag.
2. Fit on smallest part of the neck with two fingers of slack.
3. Secure the end of the collar in the belt-keeper



Re-using an active tag on a new animal

Assign the active tag to a new animal in the usual way. It will automatically be unassigned from the old animal.

Putting an active tag to sleep to conserve the battery

Hold the magnet over the magnet icon until the LED flashes ONCE.



Marking an animal as 'sold' or 'died'

1. Log in to Datamars Livestock monitoring.livestock.datamars.com
2. In the **Animals** screen, select an animal, click ... then **Mark as sold** or **died**

